

Sharpsguard® eco yellow 1

Product Code: DD477YLR

Our Sharpsguard® eco yellow range is used for the disposal of sharps that have been contaminated with medicines. Eco yellow 1 is ideal for transporting to the bedside in a POUDES® tray. Its push-down, one-piece lid incorporates a needle remover.

The container is:

- manufactured using 40% recycled material
- certified to ISO23907:2019 and UN3291
- compliant with HTM 07-01 guidance for the safe management of healthcare waste by having a yellow base and a yellow lid
- Syringe Sizes: 1, 2, 5, 10ml
- No. of Syringes*: Up to 20
- Assembly: Pre-assembled
- Accessories available:
 - DD487 POUDES® small wall bracket
 - DD497 POUDES® round leg trolley bracket
 - DD501 POUDES® square leg trolley bracket
 - DD505 POUDES® medi rail bracket
 - PT199 POUDES® community tray
 - PT199D POUDES® combination tray



Usable Capacity (l)	Product Weight (g)	Product Dimensions Assembled (max excluding bracket handles)			Aperture Dimensions (max)		
		Height (mm)	Width (mm)	Diameter / Depth (mm)	Width (mm)	Depth (mm)	Handle
0.84	154	119	204	68	70	32	Yes

Carton Dimensions					
Carton Quantity	Height (mm)	Width (mm)	Depth (mm)	Carton Weight Empty (g)	Carton Weight & Containers (kgs)
30	365	429	314	410	5.04

*based on 1ml Insulin needle with a 29g needle attached



For more information, contact 0333 060 2836 or visit stericycle.co.uk.

© 2022 Stericycle, Inc. All rights reserved.

