

Sharpsguard® eco yellow 11.5

Product Code: DD476YLR

Our Sharpsguard® eco yellow range is used for the disposal of sharps that have been contaminated with medicines. Eco yellow 11.5 is a multi-purpose container featuring a robust pail handle for easy transportation. A wider aperture allows easy insertion of larger items.

The container is:

- manufactured using 45% recycled material
- certified to ISO23907:2019 and UN3291
- compliant with HTM 07-01 guidance for the safe management of healthcare waste by having a yellow base and a yellow lid
- **Syringe Sizes:** 1, 2, 5, 10, 20, 50ml
- **No. of Syringes*:** Up to 560
- **Assembly:** Lid and base supplied separately
- **Accessories available:**
 - DD488 POUDS® large wall bracket
 - DD465 POUDS® mobile holder
 - DD465K POUDS® kickabout 11.5



Usable Capacity (l)	Product Weight (g)	Product Dimensions Assembled (max excluding bracket handles)			Aperture Dimensions (max)		
		Height (mm)	Width (mm)	Diameter / Depth (mm)	Width (mm)	Depth (mm)	Handle
9.27	622	249	N/A	275	148	48	Yes

Carton Dimensions					
Carton Quantity	Height (mm)	Width (mm)	Depth (mm)	Carton Weight Empty (g)	Carton Weight & Containers (kgs)
20	566	943	283	800	13.23

*based on 1ml Insulin needle with a 29g needle attached



For more information, contact 0333 060 2836 or visit stericycle.co.uk.

© 2022 Stericycle, Inc. All rights reserved.

