

Sharpsguard® eco yellow 2.5

Product Code: DD472YLR

Our Sharpsguard® eco yellow range is used for the disposal of sharps that have been contaminated with medicines. Eco yellow 2.5 is a general-purpose sharps container, providing maximum capacity in a minimum working space.

The container is:

- manufactured using 44% recycled material
- certified to ISO23907:2019 and UN3291
- compliant with HTM 07-01 guidance for the safe management of healthcare waste by having a yellow base and a yellow lid
- **Syringe Sizes:** 1, 2, 5, 10, 20ml
- **No. of Syringes*:** Up to 120
- **Assembly:** Lid and base supplied separately
- **Accessories available:**
 - DD487 POUDES® small wall bracket
 - DD497 POUDES® round leg trolley bracket
 - DD501 POUDES® square leg trolley bracket
 - DD505 POUDES® medi rail bracket
 - PT199D POUDES® combination tray
 - PT299 POUDES® tray 2.5



Usable Capacity (l)	Product Weight (g)	Product Dimensions Assembled (max excluding bracket handles)			Aperture Dimensions (max)		Handle
		Height (mm)	Width (mm)	Diameter / Depth (mm)	Width (mm)	Depth (mm)	
2.21	232	166	N/A	176	77	46	Yes

Carton Dimensions					
Carton Quantity	Height (mm)	Width (mm)	Depth (mm)	Carton Weight Empty (g)	Carton Weight & Containers (kgs)
48	346	663	343	550	11.68

*based on 1ml Insulin needle with a 29g needle attached



For more information, contact 0333 060 2836 or visit stericycle.co.uk.

© 2022 Stericycle, Inc. All rights reserved.

