

Sharpsguard® eco yellow 22

Product Code: DD475YLR

Our Sharpsguard® eco yellow range is used for the disposal of sharps that have been contaminated with medicines. Eco yellow 22 is ideal for theatres, laboratories and sluice rooms with high volumes of hazardous waste. The wide aperture allows economical and convenient disposal of large items.

The container is:

- manufactured using 43% recycled material
- certified to ISO23907:2019 and UN3291
- compliant with HTM 07-01 guidance for the safe management of healthcare waste by having a yellow base and a yellow lid
- **Syringe Sizes:** 1, 2, 5, 10, 20, 50ml
- **No. of Syringes*:** Up to 1050
- **Assembly:** Lid and base supplied separately
- **Accessories available:**
 - DD488 POUDS® large wall bracket
 - DD465 POUDS® mobile holder
 - DD465K POUDS® kickabout 22



Usable Capacity (l)	Product Weight (g)	Product Dimensions Assembled (max excluding bracket handles)			Aperture Dimensions (max)		Handle
		Height (mm)	Width (mm)	Diameter / Depth (mm)	Width (mm)	Depth (mm)	
17.86	943	343	N/A	327	148	58	Yes

Carton Dimensions					
Carton Quantity	Height (mm)	Width (mm)	Depth (mm)	Carton Weight Empty (g)	Carton Weight & Containers (kgs)
10	326	1148	323	780	10.21

*based on 1ml Insulin needle with a 29g needle attached



For more information, contact 0333 060 2836 or visit stericycle.co.uk.

© 2022 Stericycle, Inc. All rights reserved.

